

MGT Pump

Pumped Multi-Gas Detector

Up to six gas measurements
in a single instrument



The MGT Pump is capable of detecting up to six gases simultaneously. The internal pump activates on start up and draws the gas sample across the sensors.



Easy to use

- Self test and self monitoring
- Easy check & calibration
- 30EA of Bump tests/ Events/ Cal.
- Configuration via exclusive IR Link

Flexible

- Simultaneously detect 6 different gases
- 30 hours use from one 4-6 hour charge (P-type)
- 5 days use from one 4-6 hour charge (N-type)

The MGT Pump has four sensor slots capable of accepting a combustible gas sensor such as a catalytic pellistor or NDIR for methane. The NDIR sensor can also be a dual methane and carbon dioxide sensor. Combining the dual carbon monoxide hydrogen sulphide with another toxic gas sensor provides another three gas measurements. The additional toxic gas sensor can be either electrochemical or PID. In addition to these there is a slot for an oxygen sensor to detect low and high levels of oxygen.



Standard Package Contents

- Detector complete with stainless steel alligator clip
- Adaptor with 5 plugs
- Calibration cap
- Quick disconnect
- Dust filters (5 pieces)
- PTFE Filters (5 pieces)
- User manual
- Calibrate test certification

Model	MGT Pump
Self-test	Full function self-test upon activation Sensors, battery, Pump and circuitry – continuous
Detecting Method	Pump
Display	LCD Display (back-light)
Alarm	LED flashing, Sound, Vibration
Sensors	Electrochemical Combustible(LEL): IR or Pellistor, PID
IP	67
Temperature	-20°C ~ +50°C
Humidity	5~95% RH(Non-condensing)
Battery	Rechargeable Li-ion power supply unit (2000mAh)
Charge Time	4-6 hours
Battery Life	IR: 5 days continuous / Pellistor: 30 hours continuous
Low Battery Alarm	20 minutes / 10 minutes / 5 minutes expired: 5 long beeps/flashes and then Off is displayed
Warranty	2 years
Weight	400g
Size	77(w) × 147(h) × 54(d)mm
Optional Accessory	MGT-IR Link, Docking Station
Approval	ATEX Ex da ia IIC T4 Ga; Ex da ia IIB T4 Ga Ex ia IIC T4 Ga; Ex ia IIB T4 Ga Ex db ia IIC T4 Ga; Ex db ia IIB T4 Ga *depending upon configuration

Configuration options

Version	O2	CO/H2S	CH4 (Catalytic)	CH4 IR	CH4/CO2 IR	PID	SO2	CL2	NH3	H2	NO2
MGTP001	X	X	X								
MGTP002	X	X	X			X					
MGTP003	X	X	X				X				
MGTP004	X	X	X					X			
MGTP005	X	X	X						X		
MGTP006	X	X	X							X	
MGTP007	X	X									X
MGTP008	X	X		X							
MGTP009	X	X		X		X					
MGTP010	X	X		X			X				
MGTP011	X	X		X				X			
MGTP012	X	X		X					X		
MGTP013	X	X		X						X	
MGTP014	X	X		X							X
MGTP015	X	X			X						
MGTP016	X	X			X	X					
MGTP017	X	X			X		X				
MGTP018	X	X			X						X
MGTP019									X		
MGTP020	X	X			X					X	
MGTP021	X	X			X			X			